TONOSAFE

Instructions for use

Look closer. See further.

Tonosafe® Sterile Disposable Prisms with Tonosafe® Holder

 These devices are only to be used by healthcare professionals such as ophthalmologists, optometrists, ophthalmic technicians, nurses, or emergency physicians who have been trained in Tonometry procedures and who have read and understood these instructions for use.

INSTRUCTIONS FOR USE:

Contents:

			
5805100	x100	x5*	x1
5805091	x20	x1*	x1

* Sterile for first use only.

Indications for use:

Indicated to prevent cross infection between the eyes of different patients requiring applanation tonometry.

Intended purpose:

Tonosafe®:

A non-active, single-use accessory intended for use by healthcare professionals such as ophthalmologists, optometrists, ophthalmic technicians, nurses, or emergency physicians in either a clinical, emergency, ambulatory or domestic setting. It is intended to be connected to the Tonosafe® Holder which in turn is connected to a Perkins or Goldmann type Tonometer. It is used to prevent possible microbial cross-contamination between patients and help transmit the intraocular pressure from the eye to the Tonometer, which gives a measurement that helps in the determination of intraocular hypertension. Tonosafe® is supplied sterile ready for use and is intended to be non-contact with the user. It has direct patient contact to the mucosal membrane of the eyes for up to three applications during one procedure. The total patient exposure is determined as having a cumulative transient patient contact, i.e., less than one minute. It is recommended to wear PPE or follow facility guidelines when carrying out this procedure to prevent accidental contact and cross contamination. The Tonosafe® can be used on any patient population and is only contraindicated where medical assessment by the healthcare professional has deemed it unsafe in the case of infection, trauma or those with special needs. It is to be handled with care and stored safely to prevent damage from excessive force or atmospheric conditions.

Tonosafe® Holder:

A non-active, limited reuse accessory intended for use by healthcare professionals and designed to be connected to the Tonosafe® at the distal end and a Perkins or Goldmann type Tonometer at the proximal end. The Tonosafe® holder can be used for up to 20 procedures to secure the Tonosafe® to the Tonometer. It is supplied either sterile or clean and ready for use. It is intended to be non-contact with the patient and user, it is recommended to wear PPE or follow facility guidelines when carrying out this procedure to prevent accidental contact and cross contamination. The Tonosafe® Holder can be used on any patient population that has been risk assessed and deemed suitable for use with the Tonometer and is only contraindicated where medical assessment by the healthcare professional has deemed it unsafe in the case of infection, trauma or those with special needs. It is to be handled with care and stored safely to prevent damage from excessive forces or atmospheric conditions.

Side effects


See Tonometer IFU for side effects on using the Tonometer. The use of the Tonosafe may result in the following side effects:

- Anxiety before and during the procedure
- Pain, irritation, itching, tearing or other discomfort to the eye
- Blurred or discoloured vision from anaesthetic drops or fluorescein dye or tearing
- Bruise/contusion, corneal scarring, abrasion or other eye injury
- Infection and inflammation
- If procedure lasts too long drying can occur on the corneal epithelium of both eyes
- Topical anaesthetic agent and/or fluorescein dye may cause stinging effect on the eye
- Possible allergic reaction to anaesthetic drops or fluorescein.

Warnings, risks and other information for safety:

- Read these instructions carefully before first use. Failure to comply with these instructions may result in damage to the device and/or pose a risk to the patients and/or users.
- These instructions must be read in conjunction with the applicable Tonometer instructions for use.
- Tonosafe® devices are only for use as described in these instructions for use.
- These devices are intended to be used in a clinical, emergency, ambulatory or domestic setting.
- These devices are only to be used by healthcare professionals such as ophthalmologists, optometrists, ophthalmic technicians, nurses, or emergency physicians who have been trained in Tonometry procedures and who have read and understood these instructions for use.
- The devices are only to be handled and stored as described in these instructions, for use by facility approved and trained personnel.
- It is recommended to wear PPE or follow facility guidelines when carrying out this procedure. Training in the instructions for use is the responsibility of the clinical institute and the user.
- The Tonosafe® are for single use only.
- DO NOT touch the plane surface of the Tonosafe®.
- Tonosafe® re-use and/or re-sterilisation (including autoclaves) can result in structural damage to the device or chemical residues causing physical harm to the patient's eye and/or the risk of an incorrect prognosis. Infectious diseases can be transmitted through reuse, among which are infections that can be introduced through the mucous membranes of the eye (conjunctiva). These include viruses and bacteria that can cause conjunctivitis (e.g., adenovirus, herpes simplex, Staphylococcus aureus) and viruses that can cause systemic infections, including bloodborne viruses (e.g. hepatitis B and C viruses, human immunodeficiency virus, herpes viruses, and rhinoviruses).
- DO NOT clean prism tip with alcohol wipes or any other solution as this will cause damage to the device. The Tonosafe® are supplied sterile and ready for use.
- Incorrect use may result in scratches and/or other harm to the cornea, and/or incorrect readings.
- These devices are not UKCA or CE marked for veterinary use.
- Any serious incident that has occurred in relation to the device must be reported to the manufacturer and the Regulatory/Competent Authority of the country in which the user and/or patient is established.
- These devices do not contain latex.

Important:

-  Do not reuse
-  Do not reprocess the Tonosafe®
-  Read the instructions carefully before first use
-  Federal law restricts this device to sale by or on the order of a physician or practitioner

Handling and storage:

Store in a clean, dry place at all times to maintain safety and performance characteristics. If the packaging is damaged, soiled or unintentionally opened before use, do not use. If received damaged contact Haag-Streit UK.

-  These are delicate medical devices and must be handled with care at all times.
-  Do not use if the packaging has been damaged. If received damaged contact Haag-Streit UK.
-  Keep Dry
-  Store out of direct sunlight.

- ⚠ Before use:**
- Ensure that no damage has occurred after manufacturing and during any handling and storage by checking for signs of any physical damage, surface irregularities, imperfections or contaminants. Contact Haag-Streit UK in the event any Tonosafe® or Tonosafe® Holder has such issues, these devices should not be used or reprocessed but replaced.
 - It is recommended to wear PPE or follow facility guidelines when carrying out this procedure.
 - Avoid contact with the applanation surface to prevent cross contamination.
 - Ensure the assembly of the Tonometer, Tonosafe® Holder and Tonosafe® is secure before use.
 - Ensure patient is fully relaxed prior to taking a measurement.
 - Dry off any excess tearing from the eye to prevent getting a bad measurement.
 - If the patient squints or blinks too much, use a cotton bud or fingers to help keep the eyelid open. Make sure not to apply pressure to the globe as this will give a false reading.

⚠ Contraindications:
 These devices are contraindicated where medical assessment by the healthcare professional has deemed it unsafe in the case of infection, disease or trauma. For example: surgical procedures carried out on, or other injuries to the cornea, ruptured globe, unhealed corneal abrasions or ulcers, corneal scars or elevated eyes with astigmatism. They are also contraindicated where there is an inability to administer, or the patient is allergic to the topical anaesthetic, or otherwise when medical assessment deems it unsafe in the case of patients with special needs.

For contraindications for use of the Tonometer, refer to Tonometer instruction manual.

⚠ Cleaning and maintenance:

- Tonosafe® are single use devices and DO NOT require cleaning or maintenance.
- Tonosafe® Holders must be decontaminated between use as described below, discard after 20 uses.
- If the Tonosafe® Holder is lost or becomes damaged, contact Haag-Streit UK for replacements.

Tonosafe® Holders should be decontaminated between use by following the instructions below compiled in accordance with ISO 17664.

In the event a Tonosafe® Holder is damaged or lost, contact Haag-Streit UK for a replacement.

The instructions provided below have been validated by the manufacturer of the medical device as being capable of preparing a medical device for reuse. It remains the responsibility of the processor to ensure that the processing, as actually performed using equipment, materials and personnel in the processing facility, achieves the desired result. This requires verification and/or validation and routine monitoring of the process.

Manufacturer: Haag-Streit UK Limited. Device: Tonosafe® Holders ONLY.	
WARNINGS.	The user must follow these instructions, including a. and b. in "Before Use" above, to ensure suitable cleaning and disinfection for reuse. It is recommended to wear PPE or follow facility guidelines when carrying out this procedure.
Limitations on processing	The holder must be discarded after 20 uses. Note: Each tray contains 20 Tonosafe®.
INSTRUCTIONS.	
Preparation before cleaning	If applicable, remove from Tonometer and disassemble Tonosafe® from device
Initial treatment at the point of use	Wipe using a UKCA/CE marked disinfection wipe for non-invasive medical devices and air dry.
Cleaning: Automated	Automated cleaning is not recommended.
Cleaning: Manual	As soon as possible after use, rinse the Tonosafe® Holder fully in water ensuring the lumen is also flushed, then wipe using a UKCA/CE marked disinfection wipe for non-invasive medical devices and air dry.
Disinfection	Wipe the Tonosafe® Holder with a UKCA/CE marked disinfection wipe or soak in 70% isopropyl alcohol following their manufacturer's instructions, then air dry.
Drying	Allow Tonosafe® Holder to fully air-dry before reassembly or storage.
Maintenance, Inspection and Testing	Ensure Tonosafe® Holder has not been damaged as a result of the cleaning process before reassembly. See a. and b. in "Before Use" above.
Packaging	After use, place Tonosafe® Holder back in storage tray, see Figure 3 below.
Sterilization	The Tonosafe® Holder does not require sterilisation.
Storage	See "Handling and Storage" above.
Usage	Follow local protocols for logging the number of times used.

Other Useful Symbols:

 Swiss authorised representative	 Distributor	 Importer
 EU authorised representative	 Quantity	 Unique device identifier
 Name of manufacturer	 Catalogue number	 Classed as a Medical device
 Manufactured in Great Britain	 Tonosafe (Prism)	 Tonosafe (Holder)
 Sterilised by Ethylene Oxide Single Sterile barrier system with protective packaging outside		 Use-By Date

⚠ Technical Specification:
 The Tonosafe® Sterile Disposable Prisms are designed to meet the International Standard ISO 8612:2009.

IOP range	0 mmHg to 80 mmHg (Tolerance + 5.0 mmHg)
Flatness*	±5 µm
*Over a central area of not less than 4 mm diameter, the front surface shall be a plane structure with a "peak to valley" deviation from a plane surface of less than 3.0 µm.	

Tonosafe® assembly:

Note: Inspect the Tonosafe® and Tonosafe® Holder for signs of any nonconformances after removal from the tray (see a. and b. in "Before Use" above).

1. Carefully peel back the pre-cut top strip on the lid and extract the prism holder.
2. Carefully peel back only one of the pre-cut squares to expose a single prism tip.
3. Place the Tonosafe® Holder (B) over the exposed end of the Tonosafe® (a), see Figure 1.
4. Rotate Tonosafe® Holder (B) until the slot (D) aligns with the indexing key (C) of the Tonosafe® (see Figure 2) to create a Tonosafe® Assembly.
5. Remove the Tonosafe® Assembly from tray ensuring the indexing key is aligned.

If manual manipulation of the Tonosafe® is required to ensure a secure connection, restrict contact only to the outer flange at the distal end of the Tonosafe®. Do not touch the area of appplanation. Ensure the Tonosafe® is securely attached.

6. Attach the Tonosafe® Holder securely to the Tonometer tip, ensure the flat edge of the Tonosafe® is uppermost and horizontally aligned. It is now ready for use.
7. Follow the Tonometer instructions to obtain the Intraocular Pressure (IOP) readings.
8. When tonometry is complete, remove the Tonosafe® assembly from the Tonometer. Detach the Tonosafe® by pulling it off the Tonosafe® Holder using the flange.
9. To prevent accidental contamination after use, discard the Tonosafe® (F) by following the locally approved waste procedures (Figure 4) or by putting in the Clinical Waste Bin.
10. Decontaminate the Tonosafe® Holder immediately after use by following the 'Cleaning and maintenance' section above, and then place back into the storage tray (E) (Figure 3).


Figure 1


Figure 2


Figure 3


Figure 4

Trouble Shooting:

Fault	User Action	Further Action
Unable to connect Tonosafe® to Tonosafe® Holder.	Ensure the Tonosafe® Holder is lined up properly over the centre of the Tonosafe®.	Contact manufacturer if still unable to connect.
	Check Tonosafe® and Tonosafe® Holder for damage.	Replace either or both if damaged. Contact manufacturer if still unable to connect.
Unable to connect Tonosafe® assembly to Tonometer.	Ensure the Tonosafe® Holder is lined up properly over the centre of the tip of the tonometer.	Contact manufacturer if still unable to connect.
Device surface irregularities under magnified vision.	The device is designed to be free of surface irregularities and imperfections (free from fissures, cracks and dents) when viewed with unmagnified corrected vision.	If irregularities are present under unmagnified vision, contact the manufacturer.
Inaccurate IOP measurements (low/high).	Refer to tonometer instruction manual.	Contact tonometer manufacturer.

Contact Details:

UK
Haag-Streit UK Limited,
Unit C, Woodside Industrial Estate,
Dunmow Road, Bishop's Stortford,
Hertfordshire, CM23 5RG,
United Kingdom
info@haag-streit-uk.com
www.haagstreituk.com

EC REP Haag-Streit Deutschland GmbH,
Rosengarten 10,
22880 Wedel,
Germany
Tel: +49 4103 709 03

CH REP Haag-Streit AG,
Gartenstadtstrasse 10,
3098 Koeniz,
Switzerland.
Tel: +41 31 978 01 11
sales@haag-streit.com
info@haag-streit.com
www.haag-streit.com

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UK Haag-Streit UK Limited,