
⚠ These devices are only to be used by healthcare professionals such as ophthalmologists, optometrists, ophthalmic technicians, nurses, or emergency physicians who have been trained in ophthalmology.

Read the instructions for use carefully before commencing use of this device. It contains important information regarding the safety of the user and patient.

**INSTRUCTIONS FOR USE:**

**Contents:**

Rounded end of Schirmer Tear Test strip (Patient contact)      The opposite end of the Schirmer Tear Test strips. Each pack contains one strip with an angled end one with a squared end.

**Indications for use:**

The Schirmer Tear Test (Sterile) is a small absorbent strip of filter paper, supplied sterile, to be inserted under the patient's lower eyelid to stimulate and evaluate the formation of tears.

**Intended purpose:**

The Schirmer Tear Test Strip is indicated on patients who require assessment in cases of suspected lacrimal gland dysfunctions regarding the volume or composition of tear fluid. The Schirmer Tear Test assesses the aqueous component of the ocular tear film and can help identify aqueous tear film deficiencies. They are intended to be used on patients of any gender, ethnicity, or age after medical assessment. Caution is required during use in the case of advanced age, illness, infection, ophthalmic surgery, trauma or in those patients with special needs. They come into transient contact with the patient's tears and the palpebral conjunctiva of the lower eyelid and the bulbar conjunctiva of the eye; both eyes are normally tested. They are sterile, single use devices and intended to be used in a clinical or domestic setting by a qualified ophthalmic healthcare professional trained in their use.

They can be used with anaesthesia (Schirmer II test) as indication of the basal tear secretion or without (Schirmer I test) as an indication of basal and reflex tear secretion

**Side effects:**

The use of the Schirmer Tear Test may result in the following side effects:

- Anxiety before and during the procedure
- Pain, irritation, itching, tearing or other discomfort to the eye
- Infection and inflammation
- If procedure lasts too long, drying can occur on the corneal epithelium of both eyes.

⚠ **Warnings, risks, and other information for safety:**

- Read these instructions carefully before first use. Failure to comply with these instructions may result in damage to the device and/or pose a risk to the patients and/or users.
- Schirmer Tear Test devices are only for use as described in these instructions for use.
- These devices are intended to be used in a clinical, emergency, ambulatory or domestic setting.
- These devices are only to be used by healthcare professionals such as ophthalmologists, optometrists, ophthalmic technicians, nurses, or emergency physicians who have been trained in ophthalmic procedures and who have read and understood these instructions for use.
- The devices are only to be handled and stored as described in these instructions, for use by facility approved and trained personnel.
- Personal Protective Equipment (PPE) must be worn when using these devices to prevent cross infection and microbial hazards. Follow facility approved procedures for the use of PPE.
- Training in the instructions for use is the responsibility of the clinical institute and the user.
- Schirmer Tear Test devices are for single use only. It is recommended to use one strip per eye to prevent possible microbial cross-contamination between the patient's eyes.
- DO NOT touch the rounded wick end of the Schirmer Tear Test.
- Schirmer test should not be the sole test for diagnosing dry eye disease. The test should be used in conjunction with other diagnostic procedures to confirm diagnosis.
- Re-use and/or re-sterilisation (including autoclaves) can result in structural damage to the device or chemical residues causing physical harm to the patient's eye and/or the risk of an incorrect prognosis. Infectious diseases can be transmitted through reuse, among which are infections that can be introduced through the mucous membranes of the eye (conjunctiva). These include viruses and bacteria that can cause conjunctivitis (e.g. adenovirus, herpes simplex, Staphylococcus aureus) and viruses that can cause systemic infections, including bloodborne viruses (e.g. hepatitis B and C viruses, human immunodeficiency virus), herpes viruses, and rhinoviruses.
- DO NOT clean the Schirmer Tear Tests with alcohol wipes or any other solution as this will cause damage to the device. The Schirmer Tear Test are supplied sterile and ready for use.
- Incorrect use may result in scratches and/or other harm to the cornea.
- These devices are not UKCA or CE marked for veterinary use.
- Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the Regulatory/Competent Authority of the country in which the user and/or patient is established.
- These devices do not contain latex.

⚠ **Important:**

- ⊘ Do not reuse
- ⌛ Do not reprocess the Schirmer Tear Test

📖 Read the instructions carefully before first use      **Rx Only** Federal law restricts this device to sale by or on the order of a physician or practitioner

⚠ **Handling and storage:**

Store in a clean, dry place at all times to maintain safety and performance characteristics. If the packaging is damaged, soiled or unintentionally opened before use, do not use. If received damaged contact Haag-Streit UK.

🍷 These are delicate medical devices and must be handled with care at all times.      ⚠ Do not use if the packaging has been damaged. If received damaged contact Haag-Streit UK.

☂ Keep Dry      ☀ Store out of direct sunlight.

⚠ **Before use:**

- The examination is to be performed at the medical professional's discretion and consideration should be taken in cases of active infections and lesions of the eye
- Ensure that no damage has occurred after manufacturing and during any handling and storage by checking for signs of any physical damage, imperfections, or contaminants. Contact Haag-Streit UK in the event of any such issues. These devices should not be used or reprocessed but replaced
- It is recommended to wear protective medical gloves when carrying out the procedure
- Avoid contact with the rounded wick end to prevent cross contamination
- Ensure patient is fully relaxed prior to carrying out the procedure.

**Contraindications:**

- These devices are contraindicated where:
- Medical assessment by the healthcare professional has deemed it unsafe in the case of infection, disease, or trauma. For example, surgical procedures carried out on, or other injuries to the cornea, ruptured globe, unhealed corneal abrasions or ulcers, corneal scars or elevated eyes with astigmatism
  - When medical assessment deems it unsafe in the case of patients with special needs.
  - When a patient has just received eye drops

⚠ **Cleaning and Maintenance:**

Schirmer Tear Test Strips are single use devices and **DO NOT** require cleaning or maintenance. The Haag-Streit Schirmer Tear Tests are made of Cellulose paper

Please take note of the following restrictions:

- Do not re-sterilise
- Do not expose to temperatures above 60°C

The above-mentioned instruction manual, the brief instructions and additional information can always be found on the website: [www.haagstreituk.com](http://www.haagstreituk.com).

**Environmental Conditions:**

Device is designed to be used, transported, and stored in ambient conditions.

⚠ Do not use the device if it is exposed to environmental conditions outside of those specified above. Contact Haag-Streit UK for further guidance.

**Other Useful Symbols:**

|   |   |                                |
|---|---|--------------------------------|
| CH REP Swiss authorised representative  | Distributor   | Importer                       |
| EC REP EU authorised representative   | QTY Quantity  | UDI Unique device identifier   |
| Name of manufacturer  | REF Catalogue number  | MD Classed as a Medical device |
| Manufactured in Great Britain   | Use-By Date   |                                |
| Sterilised by Gamma Radiation Single Sterile barrier system with protective packaging outside | Sterilised by Gamma Radiation Single Sterile barrier system |                                |

**Warranty & product liability**

Haag-Streit UK products must be used only for the purposes and in the manner described in the documents distributed with the product.

Improper handling can damage the product. This would void all warranty claims.

**Procedure**

The Schirmer Tear Test strips serve as an auxiliary tool to evaluate lacrimation. These are typically employed during a comprehensive clinical examination, where all environmental and clinical variables are considered. The results from the Tear Test strip are useful for the evaluation of a patient's suitability for wearing contact lenses. This device is intended for exclusive use by medically trained health care professionals.

- Inform the patient about the simplicity and necessity of this test in understanding their tear production. Clarify that while the test is harmless and painless, it could cause slight discomfort in some instances.
- Seat the patient in an examination chair, ideally with a headrest, in a room with attenuated lighting.

Separate the packaging pouch from the opposite end, avoiding the wick area. Carefully grasp the rounded wick end over the pouch to prevent contamination. Proceed to bend the rounded wick end of the test strip at the marked 'V' indentation (Fig 1a) to form a right angle (Fig 1b).

Remove the strip from the packaging, making certain not to touch the rounded wick end to avoid contamination.

Instruct the patient to look upward while you gently pull down the lower eyelid. Place the rounded, bent end of the sterile strip on the lower fornix so that the bend overlaps the lid margin.

Take precautions to avoid causing irritation while positioning the test strip, as this could trigger excessive tear flow and compromise the test results. If necessary, topical anaesthetics may be utilised to prevent irritation.

After a period of 5 minutes, assess the length of the moistened section of the test strip (Fig 1c). Results will differ dependent on whether Schirmer I or Schirmer II is used.

**Guidelines to be mindful of during standard test routine:**

- The test should be carried out prior to any other intervention involving the eye
- Ensure that the strip does not come into contact with the cornea, as this would invalidate the test results
- The patient may blink normally during the test, although some individuals may prefer to keep their eyes closed.

**Disposal**

To prevent accidental contamination after use, discard the Schirmer Tear Tests by following locally approved waste procedures or by putting in the clinical waste bin.

**Trouble Shooting:**

| Fault  | User Action   | Further Action   |
|--|---|--|
| Device surface irregularities under magnified vision | The device is designed to be free of surface irregularities and imperfections when viewed with unmagnified corrected vision | If irregularities are present under unmagnified vision, contact the manufacturer |

**Contact:**

Haag-Streit UK Limited,  
Unit C, Woodside Industrial Estate,  
Dunmow Road, Bishop's Stortford,  
Hertfordshire, CM23 5RG,  
United Kingdom  
info.unitedkingdom@haag-streit.com  
www.haagstreituk.com

**EC REP** Haag-Streit Deutschland GmbH,  
Rosengarten 10,  
22880 Wedel,  
Germany  
Tel: +49 4103 709 02

**CH REP** Haag-Streit AG,  
Gartenstadtstrasse 10,  
3098 Koeniz,  
Switzerland.  
Tel: +41 31 978 01 11  
sales@haag-streit.com  
info@haag-streit.com  
www.haag-streit.com

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