




UPGRADE KITS

Extend the lifetime of your device

Look closer. See further.





Haag-Streit Upgrade Kits — Extend the lifetime of your device

Haag-Streit is world-renowned for developing, manufacturing, and distributing high-quality products and solutions backed by outstanding customer service.

To prolong the lifetime of your device – and increase your return on investment – we are pleased to offer you our innovative range of Upgrade Kits.

What is a Haag-Streit Upgrade Kit?

An Upgrade Kit is a convenient package containing compatible components that can enhance specific Haag-Streit devices. Your choice of Upgrade Kit will depend on your particular processes and activities – some will enable you to implement new techniques and diagnostic assessments, while others are designed to make your daily activities even more comfortable. Regardless of which Upgrade Kit you choose, you can be assured that the upgrade will be implemented by Haag-Streit Service Engineers to ensure our “gold-standard” service.



Benefits of a Haag-Streit Upgrade Kit

On-site or Service Center – we take care of it

We appreciate that time is essential to you and, therefore, make every effort to keep downtime to a minimum. Most upgrades can be carried out swiftly and conveniently at your site.

On occasions when the device is required to travel to our Service Center, we are pleased to offer loaner devices, while the upgrade is carried out, so that you are never without a Haag-Streit device.

Regardless of the Upgrade Kit you choose, relax knowing that your upgrade is handled by a Haag-Streit Accredited Service Engineer who is fully trained to make the necessary adjustments.

Upgrading: the smart, sustainable choice

By choosing to upgrade, you automatically extend the lifetime of your device, decreasing the need for a complete product replacement, which is costly to the planet – and your budget. In addition, some Haag-Streit Upgrade Kits inherently reduce negative environmental impacts.

Read on to learn about the different Upgrade Kits available.

Lenstar Upgrade Kits

Benefit from enhanced illumination plus other value-adding upgrades

Upgrade your Lenstar from “green light” to “white light” and enjoy all the benefits of reliable, high-quality color images when taking white-to-white measurements. This illumination upgrade provides greater detail of anatomical landmarks and orientation in the surgical planner and in the OR.

Moreover, this Upgrade Kit provides a gateway into additional enhancements and accessories, such as the ability to add a T-Cone. The T-Cone topography accessory supports the identification of irregular astigmatism and allows double-checking on the axis position. Together with enhanced illumination, this upgrade enables well-informed IOL selection and increased safety of toric IOL surgery.

Plus, when you upgrade, our Accredited Service Engineers carry out a complete check of the measurement units in your device – returning a factory-calibrated device that’s as good as new.

Why upgrade?

- Benefit from high resolution and color-accurate pictures
- Take advantage of additional Lenstar upgrades & optional accessories (T-Cone, APS)
- Fitting: Service Center.

Ask us for:

- Lenstar with “white light” (part number: 1024646S)
- Lenstar with “white light” plus T-cone (part number: 1024647S).





Lenstar APS – Axial length measurement made simple

Transform your biometry examinations with our Automated Positioning System (APS) for the Lenstar. This Upgrade Kit motorizes the Lenstar's measurement process for effortless maneuvering and focus. As soon as the patient looks into the Lenstar, the APS tracks the patient's eye, and the measurement is taken with one click of the joystick. Taking quick, accurate measurements enables easy delegation and efficient patient flow in a busy practice.

Please note that you must firstly upgrade your Lenstar illumination from "green light" to "white light" to take advantage of this Upgrade Kit.

Why upgrade?

- Motorizes the Lenstar measurement process
- Speeds up biometry examinations
- Ideal for examining younger patients as part of myopia management
- Fitting: On-site or Service Center.

Ask us for:

- Upgrade Kit Lenstar 900 APS (part number: 1024648S).

Slit lamp Upgrade Kits

Slit lamp LED – the cleaner, greener, cost-saving solution

Our slit lamp LED Upgrade Kits are very popular with customers, and for good reason: switching older tungsten bulbs to LED can save you money, as well as reducing negative impacts on the environment. Long-lasting LEDs provide stable, brilliant illumination for up to 150 times longer than tungsten – reducing the time and cost of routine bulb changes, while also providing an energy saving of around 60%.

Available for the BQ 900, BP 900 and BM 900, simply tell us your slit lamp model and we will talk you through the configurations available – considering your preferred power supply, illumination settings, and headrest type.

Moreover, this Upgrade Kit also automatically extends the slit height from 8mm to 14mm, making it possible to carry out specific diagnostic procedures, such as those required for laser surgery.

Why upgrade?


- Improved illumination, allowing finer detail for clinical diagnosis
- Eliminates the need for routine bulb changes, saving you money
- We configure your LED Upgrade Kit to suit your needs
- Fitting: On-site.

Ask us for:

- LED lighting for BQ 900 – black (part number: 1024649S)
- LED lighting for BQ 900 – light grey (part number: 1024650S)
- LED lighting for BP 900 – black (part number: 1024651S)
- LED lighting for BM 900* – black (part number: 1020884S).

* LED upgrade only - does not include 14mm extension





Slit lamp joystick – enjoy the same ease of use for even longer

Whether you are in your early career or a seasoned ophthalmologist: your slit lamp is an essential tool that undergoes intense manual handling over its lifetime. Switch out your rubber joystick with a metal replacement and enjoy the same focus and ease of use with additional durability. Very quick to replace, the new joystick isn't just easy to clean and disinfect – it looks great, too.

Why upgrade?

- Sleek, long-lasting metal finish
- Easy to remove
- Easy to clean & compliant with hospital disinfection procedures
- Fitting: On-site.

Ask us for:

- Slit lamp joystick (part number: 1024653S).

[Find out more...](#)

Interested to learn more about how upgrading can enhance the lifespan of your Haag-Streit device? To discuss any of the Upgrade Kits shown in this brochure, please contact your Support Agent who will be delighted to help.

HAAG-STREIT UK

Unit C, Woodside Industrial Estate,
Dunmow Road,
Bishop's Stortford,
CM23 5RG
Phone +44 (0)1279 414969
info.unitedkingdom@haag-streit.com
www.haagstreituk.com

